

C A B I N E T P A N E L S



SELECTIONS

Since cabinetry is often the dominant design element in a room, the color you select will set the overall tone and style of your dental practice. And, because of the demands on equipment in the healthcare environment, it is good to know that the materials being used to create the finish on your cabinets are as durable as they are beautiful.

Integra® Dental Casework

Membrane-pressed thermofoils, also known as 3D laminates, seal the Integra® MDF cabinet panels. This dimensionally flexible material provides a continuous, seamless surface, eliminating the need for edgebanding or special edge treatments while offering a high level of scratch and stain resistance.

Midmark offers the Integra® line of cabinetry in 13 beautiful, standard colors: 8 woodgrains, 2 solids and 3 metallics.

Artizan® Dental Furniture

Thermally-fused melamine and high pressure laminate (HPL) are used to cover Artizan® cabinet panels. These materials are applied at high temperatures and pressure to form a durable finish. Matching 3mm edgebanding is applied to completely seal the corners from moisture.

Midmark offers the Artizan® line of cabinetry in 7 beautiful, standard colors: 5 woodgrains and 2 solids. Custom colors are available by request.

To assist in your selection of Midmark countertops, upholstery and cabinetry panels, please access our Dental Color Selector at midmark.com.

The color swatches shown are the best representation of the original material. Actual color may vary slightly.


Integra® Colors



I - Fusion Maple*




I - Summer Apple



I - Shimmer (metallic)



A - Fusion Maple*




A - Frosty White



I - Golden Anigre




I - Crossfire Java*



I - Bark (metallic)



A - Ginger Root Maple




A - Pebble Grey*



I - Honey Maple



I - Dark Walnut




I - Shark Grey (metallic)



A - Sunset Cherry



I - Wild Cherry



I - Harmony White




A - Violin Figured Annigre*



I - Violin Figured Annigre*



I - Pebble Grey*



A - Crossfire Java*

I - Integra® Colors
A - Artizan® Colors

* These colors are available for both the Integra® and Artizan® lines. They are close matches; however, because they are different materials, they may not appear exactly the same. If you plan to purchase cabinets from both lines and want to match colors, we strongly recommend contacting Midmark Customer Service at 1-800-MIDMARK to request samples before placing an order.

We have tried our very best to accurately represent product colors in this brochure. However, due to the nature of printed media, the color of actual products may differ slightly from how they appear on your particular printed piece. In order to confirm your color selection, we strongly recommend that you contact Midmark Customer Service at 1-800-MIDMARK to request a sample before placing your order.

All trademarks and registered trademarks are the property of their respective owners.

CARB 93120.2 Phase 1 Compliant

Midmark is an ISO 13485 and ISO 9001 Certified Company.
For more information or a demonstration, contact your Midmark dealer or call: 1-800-MIDMARK Fax: 1-877-725-6495
Outside the U.S.A. call: 1-937-526-3662 Fax: 1-937-526-8214
or visit our website at midmark.com

© 2009 Midmark Corporation
60 Vista Drive P.O. Box 286 Versailles, Ohio 45380-0286
Products subject to improvement changes without notice
Litho in U.S.A. 007-0752-00 Rev. B1 (2/11)



Because we care.